www.myschlick.com

Living for Solutions.

Das Anti-Bearding-Konzept

The anti-bearding concept

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Innovation and tradition as stimulus for sustainability and quality

Comprehensive expertise

In the field of pressure and multiple-substance nozzles, the name SCHLICK is synonymous with quality and precision around the world. Yet for some time now the company has been doing much more than simply developing and manufacturing high-precision nozzles. With their extensive knowledge, SCHLICK specialists are increasingly becoming a firm fixture in the planning stage for new developments, such as production plants and manufacturing processes.

Responsible design

Trust, reliability and sustainability are amongst the fundamental company values at SCHLICK, and have been for more than 111 years. This involves understanding that design and technology must work together in harmony. The preservation of resources and environmental protection are key to the conduct of any business.

Proactive customer focus

Anyone wanting to set technical standards needs to keep their ear to the ground on the market. In addition to countless customer solutions, this approach has given rise to many trendsetting new 'Made-by-Schlick' developments. This innovative prowess is clear to see in developments such as the company's patented ABC technology® and Multispray technology.

„Schlick solutions – not designed by machines,
but created by our engineers.“

Was bedeutet Anti-Bearding?

What's anti-bearding supposed to mean?

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A comparison of the flow patterns reveals that the dust particle flow – simulated here in the form of fog – stays away from the anti-bearding cap much better than it does from the conventional flat spray cap. The ABC cap thus counteracts the conventional bearding effect.

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The air cap of two-substance nozzle with horns shows the typical product build-up (bearding). This clogs the air ducts and impairs the spraying pattern. The SCHLICK ABC cap, however, does not experience any product build-up at the nozzle orifices and thus allows for an optimal spraying pattern.

Redefining dual flat-spray nozzles

During coating, more and more product particles gradually start to accumulate on a conventional flat-spray cap due to the quantity of dust that is circulating. This not only influences the spray behaviour, but leads to the nozzle becoming clogged. This is also referred to as „beard formation“, or bearding.

The SCHLICK Anti-Bearding Cap (ABC) has been developed to counteract this bearding. Turbulence in the region of the air cap and the nozzle opening is minimised, meaning that interruptions to the process – due to the nozzle needing to be cleaned – can be almost completely prevented.

To make it possible to upscale/downscale across the entire production process, the SCHLICK ABC range has been developed for various flow rates, and covers drum sizes from 8.5“ to 24“ and above.

The spray has constant and homogeneous liquid distribution with very fine, uniform and reproducible droplet size distribution at the same time, and is characterised by a very uniform ellipse, which is achieved by the balance between the pressures – the forces of the atomising air and the pattern air.

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The optimal pressure for atomising and pattern air is controlled by the parameters quantity of liquid, density and viscosity as well as solid content: You must accommodate the respective application.



Flüssigkeitsmenge in g/min	Bohrung in mm (Viskosität < 100 mPas)	Bohrung in mm (Viskosität < 100 mPas)	Min. empfohlener Zerstäubungsluftdruck (AA) in bar (ü)
Quantity of liquid in g/min	Bore diameter in mm (viscosity < 100 mPas)	Bore diameter in mm (viscosity < 100 mPas)	Minimum recommended atomising air (AA) pressure in bar (g)
30 - 60	0.8 - 1.0	1.0 - 1.2	0.7
60 - 120	1.0 - 1.2	1.2 - 1.5	1.0
120 - 150	1.2 - 1.5	1.5 - 1.8	1.5
150 - 180	1.5 - 1.8	1.8 - 2.2	2.0

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
The bore diameter of the nozzle should be aligned dependent on the viscosity of the liquid.

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A **round jet** is produced without using pattern air. This jet is characterized by over-humidification in the centre and fine droplets in the area at the edges.

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An **extreme flat spray** is produced by an excessively high quantity of pattern air. The high volume flow density and overlapping spraying patterns lead to over-humidification and large droplets at the edges. Fine coverage in the centre of the spray is also typical.

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An **ideal ABC spray** ensures uniform droplets and is the result of the balance between atomising air and pattern air.

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The design of the SCHLICK Anti-Bearing Cap is patented and therefore legally protected.

GMP: Good Manufacturing Practise Design

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GMP: Good Manufacturing Practise Design

The SCHLICK ABC Technology® provides perfect upscaling of the production parameters from a tabletop device to laboratory systems and production.

In addition, a production nozzle has been developed both for systems with single/separate attachments for controlling pattern air (PA) and atomising air (AA), to provide all production systems with the advantages of the SCHLICK ABC Technology®.

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- ✓ *Sophisticated construction for very simple handling*
- ✓ *Simple maintenance and cleaning thanks to manual installation/removal*
- ✓ *Only seven or eight separate parts and four O-rings (depending on model)*
- ✓ *Including cleaning needle*
- ✓ *Liquid return system available*
- ✓ *Materials conform to FDA: 1.4404 (AISI 316 L), EPDM-O-rings*

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Models 930/7-1 S35 respectively S 45 and 970/7-1 S75 can be converted into a pressure nozzle with a flat-spray tip with very little effort, and then can be used for coating items with sugar solutions. The standard flat-spray tips from the 650 model range are used for this purpose.

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
Production nozzle with single/separate attachments for controlling pattern air and atomising air

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The design of the SCHLICK Anti-Bearding Cap is patented and therefore legally protected.

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
Our patented SCHLICK Anti-Bearding Cap prevents product build-up throughout the entire production process. By consistently using the ABC Technology®, this first of all allows perfect upscaling and downscaling. This is an important criterion for an optimally homogeneous and reproducible spraying result.

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Custom version of model S35 for mounting on the PCA with only two screws. All media are fed through the specially developed PCA connection. There are no further connectors; dirty corners are avoided.

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ANTI-BEARDING TO THE HIGHEST STANDARD

The Professional Coating Arm was specially developed for coating applications and is included with the patented ABC Technology® as standard.


The PCA is supplied via individually adapted attachments for control air, atomising air, pattern air and liquid (with flow and return) and functions as a supply unit for all nozzles, as no further hoses or connectors are required. This keeps dirty corners and installation work to a minimum. Fully homogeneous operation of up to eight nozzles is possible with the PCA. Its sophisticated construction drastically reduces maintenance and

set-up times. This is because no special tools are required to install or de-install the supply unit – only two screw fittings need to be loosened. Errors in nozzle spacings arising during installation are also ruled out thanks to the clearly defined dimensions of the individual blocks. The individual nozzles are fixed using two screws only. Each individual nozzle consists of just seven individual pieces plus O-rings.

The SCHLICK Professional Coating Arm is equipped with a customised attachment mechanism, enabling PCA technology to be operated in almost any system.

In order to reduce the weight, most of the components are manufactured from titanium 3.7035. The remaining parts are all made of stainless steel 1.4404 (Aisi 316 L). Standard EPDM O-rings (FDA-compliant).

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With very little effort, the PCA can be used for **coating items with sugar solutions**. The standard flat-spray tips from the 650 model range are used for this purpose.

Including separate attachments for control air, atomising air, pattern air and liquid (with flow and return). The attachment mechanism is also adapted to the relevant coater.

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PCA with model 930 7-1 S37 for use in tablet coating

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Standard solutions of the highest quality are business as usual for us, the leading nozzle manufacturer. Thanks to our high degree of vertical integration and comprehensive expertise, we of course can also produce custom-made and unique products for particular applications. No matter how small the batch or short the time-frame.

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PCA with model 930 7-1 S38, individual, separate supply of the nozzles with supply protection

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Spray arm with model 930 7-1 S68, especially for the top spray process in fluidised bed systems

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MEASURABLE SUCCESS

When it comes to finding the best solution, a comprehensive quality system helps enormously. SCHLICK's new Test & Research Center is a vital building block for achieving that, combining highly accurate laser optic measuring systems with the knowledge of proven experts. Modular and flexible test set-ups replicate complex spray and coating scenarios with 500 m² surfaces. This has resulted in test arrangements that correspond to real-life production conditions on the market.

Using state-of-the-art technology, the various spray parameters are measured, compared and analysed. Thanks to the company's large-scale vertical integration, the test results can be immediately implemented and retested again. This results in the optimum solutions for customers. The series production of new SCHLICK products and special constructions is subject to the same procedures. Take the NANO-PCA, for example: This lightweight type of PCA provides long-lasting support to the user throughout the workflow, without performance losses. The key to success lies in the detail.

Schlick
Atomizing technologies

Living for Solutions.

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Hutstraße 4

D-96253 Untersiemau/Coburg

Germany

Tel.: +49 (0) 95 65/94 81-0

Fax: +49 (0) 95 65/28 70

www.myschlick.com

info@myschlick.com

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Düsen-Schlick GmbH